



# Microneedle Kit



# 1

Make sure skin is clean and makeup free.  
Numb your skin lightly with 5% lidocaine cream by applying a thin translucent layer all over face.  
Leave on for 30 minutes, then remove with damp cloth.



Wash skin with choice of cleanser to remove excess cream, then wipe entire face with alcohol pads. Break open top of Hyaluronic Peptide ampule and uncap.

# 2

# 3



Uncap microneedle stamper cap and swab needles with alcohol, then unscrew full top, Pour contents of Hyaluronic Peptide ampule into stamper bottle. Screw microneedle stamper top back on and begin stamping cheeks, upper lip and chin with moderate pressure. Stamp forehead, under eyes and nose with lighter pressure until all serum has been deposited into skin. Recap stamper.



Wash hands. Apply Sheild Recovery balm with clean hands, sealing the skin and locking in moisture. Finish by applying SPF all over. Dispose of your microneedle stamper and supplies. Please do not store fluid or needle stamper. Can not be re-filled or re-used.

# 4

For personal and clinical use only. May not be resold or distributed.